- WebDia -

m.webdia.ch

Application mobile

- Tablettes tactiles
- Smartphones

- Inscription -

The image shows the registration screen of the WebDia mobile application. The header includes a home icon, the text 'Accueil', and a title bar '- Inscription -'. The form consists of several input fields: a dropdown menu for gender (currently showing 'M.'), followed by fields for 'nom', 'prénom', 'email', 'mot de passe', and 'confirmer mot de passe'. Below these fields is a CAPTCHA image with the text 'Click to change' and a corresponding input field labeled 'saisir le code'. There is a checkbox for accepting terms, with the text 'J'accepte les [conditions d'utilisation](#)'. A large blue 'Ok' button is at the bottom. The footer contains the text '- © 2013 - kmConcept -'.

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- © 2013 - kmConcept -

- Calcul Bolus -

The screenshot shows a mobile application interface for calculating insulin bolus. At the top, there is a navigation bar with a home icon and the text 'Accueil', and a title bar with '- Calcul Bolus -'. Below this, the interface contains several input fields: 'Nombre HdC' with the value '45', a meal selection dropdown menu currently showing 'Midi', 'X unité / X HdC' with values '1' and '12' separated by a slash, and 'Glycémie' with the value '8.5'. A blue button labeled 'Calcul injection' is positioned below the input fields. At the bottom of the screen, there is a footer with the text '- © 2013 - kmConcept -'.

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- Calcul Bolus -

The screenshot shows a mobile application interface with a light blue background. At the top, there is a header bar with the text "- Saisie quantité -". Below this, a white rounded rectangle contains the text "Vous devez saisir la quantité de: Pâtes cuites." followed by a text input field labeled "Quantité gr/ml". Below the input field are two buttons: "Ok" and "Annuler". At the bottom of the screen, there is a footer bar with the text "- © 2013 - kmConcept -".

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- Calcul Bolus -

The image displays two screenshots of a mobile application titled "Calcul Bolus".

The left screenshot shows the main input screen. At the top, there is a header bar with "Accueil", "- Calcul Bolus -", and "+ Ajout". Below this, the screen contains several input fields: "37.50 HdC" with a minus button, "Pâtes cuites" with a grid icon, "Ajout HdC" with a "nombre HdC" input field, "Midi" with a dropdown arrow, "X unité / X HdC" with "1" and "12" in input fields separated by a slash, and "Glycémie" with a "Glycémie (mmol/l)" input field. A large blue button labeled "Calcul injection" is at the bottom. The footer shows "- © 2013 - kmConcept -".

The right screenshot shows the same screen but with a modal dialog box overlaid. The dialog box is blue with white text that says "Injecter : 3.30" and has a white "Ok" button. The background screen shows "X unité / X HdC" with "1" and "15" in input fields, and "Glycémie" with "9" in the input field. The footer also shows "- © 2013 - kmConcept -".

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- Ajout HdC -

The screenshot shows a mobile application interface with a top navigation bar containing three tabs: 'Accueil' (with a home icon), '- Aliment -', and 'Bolus'. Below this, there are three buttons: 'Liste' (with a plus icon), 'Mon aliment' (with a checkmark icon), and 'Ajout HdC' (with a checkmark icon and highlighted in blue). The main content area features a form titled 'Quantité de HdC :'. Inside the form, there is a text input field containing the placeholder text 'nb. HdC'. Below the input field is a large 'Ok' button. At the bottom of the screen, there is a footer that reads '- © 2013 - kmConcept -'.

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- Calcul Bolus -

The screenshot shows the 'Calcul Bolus' app interface. At the top, there are three tabs: 'Accueil' (with a home icon), '- Calcul Bolus -', and '+ Ajout'. Below the tabs, there is a list of food items: '40.50 HdC' with a minus icon, 'Pâtes cuites' with a plus icon, and 'mon aliment 1' with a plus icon. Below this list, there is a field 'Ajout HdC' with the value '25'. Underneath, there is a dropdown menu currently showing 'Dîner'. Below the dropdown, there is a field 'X unité / X HdC' with the value '1' and a field '11'. Below that, there is a field 'Glycémie' with the value '12'. At the bottom, there is a large blue button labeled 'Calcul injection'. At the very bottom, there is a small text: '- © 2013 - kmConcept -'.

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The screenshot shows a mobile application interface for calculating a bolus. At the top, there is a navigation bar with a home icon and the text 'Accueil' on the left, and '- Calcul Bolus -' on the right. Below the navigation bar, there is a large blue modal dialog box with white text that reads: 'Votre glycémie est bloquée à 5, car il y a moins de 120 mn. que vous avez fait votre dernière injection !'. Below the text in the dialog is a white button with the text 'Ok'. Behind the dialog, the main interface is visible, showing a 'Glycémie' label, a red button labeled 'Débloquer', a text input field containing the number '5', and a large blue button labeled 'Calcul injection'. At the bottom of the screen, there is a small copyright notice: '- © 2013 - kmConcept -'.

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The screenshot shows a mobile application interface for calculating insulin bolus. At the top, there is a navigation bar with a home icon and the text 'Accueil', and a title bar with '- Calcul Bolus -'. The main content area has a light blue background. It contains several input fields and buttons: a 'Nombre HdC' field with the value '45' and a small up/down arrow icon; a 'Collation de 16h' field with a dropdown arrow; an 'X unité / X HdC' field with '1' and '15' separated by a slash; a 'Glycémie' field with the value '5'; a red 'Débloquer' button; and a blue 'Calcul injection' button. At the bottom, there is a footer with '- © 2013 - kmConcept -'.

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